

Product Datasheet

Human Transcription Factor Sox 1 ELISA Kit (orb1214142)

Description	Human Transcription Factor Sox 1, SOX1 ELISA Kit
Reactivity	Human
Range	2-70ng/ml
Tested Applications	ELISA
Target	SOX1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,023
Uniprot ID	000570
Solubility (25°C)	80ng/ml
Sensitivity	0.109ng/ml
Entrez	6656
Expiration Date	Please enquire.