

Product Datasheet

Human Secreted Frizzled Related Protein 1 ELISA Kit (orb1214145)

Description	Human Secreted Frizzled Related Protein 1, SFRP1 ELISA Kit
Reactivity	Human
Range	1-400ng/L
Tested Applications	ELISA
Target	SFRP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,386
Uniprot ID	Q8N474
Solubility (25°C)	480ng/L
Sensitivity	0.46ng/L
Entrez	6422
Expiration Date	Please enquire.