

Product Datasheet

Human Troponin T ELISA Kit (orb1214159)

Description	Human Troponin T, Cardiac Muscle, TNNT2 ELISA Kit
Reactivity	Human
Range	0.5-150ng/ml
Tested Applications	ELISA
Target	TNNT2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,924
Uniprot ID	P45379
Solubility (25°C)	160ng/ml
Sensitivity	0.29ng/ml
Entrez	7139
Expiration Date	Please enquire.