

Product Datasheet

Human Gastrokine 1 ELISA Kit (orb1214165)

Description	Human Gastrokine 1, GKN1 ELISA Kit
Reactivity	Human
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	GKN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,999
Uniprot ID	Q9NS71
Solubility (25°C)	96ng/ml
Sensitivity	0.16ng/ml
Entrez	56287
Expiration Date	Please enquire.