

Product Datasheet

Human Carbonic Anhydrase 9 ELISA Kit (orb1214166)

Description	Human Carbonic Anhydrase 9, CA9 ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	CA9
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	49,698
Uniprot ID	Q16790
Solubility (25°C)	40ng/ml
Sensitivity	0.021ng/ml
Entrez	768
Expiration Date	Please enquire.