

## Product Datasheet

### Human Carbonic Anhydrase 4 ELISA Kit (orb1214167)

|                            |                                           |
|----------------------------|-------------------------------------------|
| <b>Description</b>         | Human Carbonic Anhydrase 4, CA4 ELISA Kit |
| <b>Reactivity</b>          | Human                                     |
| <b>Range</b>               | 0.2-60ng/ml                               |
| <b>Tested Applications</b> | ELISA                                     |
| <b>Target</b>              | CA4                                       |
| <b>Storage</b>             | 2-8°C                                     |
| <b>Note</b>                | For research use only                     |
| <b>Application notes</b>   | Sandwich                                  |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates    |
| <b>MW</b>                  | 35,032                                    |
| <b>Uniprot ID</b>          | <b>P22748</b>                             |
| <b>Solubility (25°C)</b>   | 64ng/ml                                   |
| <b>Sensitivity</b>         | 0.093ng/ml                                |
| <b>Entrez</b>              | <b>762</b>                                |
| <b>Expiration Date</b>     | Please enquire.                           |