

Product Datasheet

Human Bone Morphogenetic Protein 3 ELISA Kit (orb1214178)

Description	Human Bone Morphogenetic Protein 3, BMP3 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	BMP3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	53,372
Uniprot ID	P12645
Solubility (25°C)	2400ng/L
Sensitivity	4.96ng/L
Entrez	651
Expiration Date	Please enquire.