

Product Datasheet

Human Bone Morphogenetic Protein 1 ELISA Kit (orb1214179)

Description	Human Bone Morphogenetic Protein 1, BMP1 ELISA Kit
Reactivity	Human
Range	5-1000ng/L
Tested Applications	ELISA
Target	BMP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	111,249
Uniprot ID	P13497
Solubility (25°C)	1280ng/L
Sensitivity	2.67ng/L
Entrez	649
Expiration Date	Please enquire.