

Product Datasheet

Human Fascin ELISA Kit (orb1214192)

Description	Human Fascin, FSCN1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	FSCN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	54,530
Uniprot ID	Q16658
Solubility (25°C)	4800ng/L
Sensitivity	11.47ng/L
Entrez	6624
Expiration Date	Please enquire.