

Product Datasheet

Human Interleukin-23 Receptor ELISA Kit (orb1214196)

Description	Human Interleukin-23 Receptor, IL-23R ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	IL23R
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	71,722
Uniprot ID	Q5VWK5
Solubility (25°C)	9600ng/L
Sensitivity	16.67ng/L
Entrez	149233
Expiration Date	Please enquire.