

Product Datasheet

Human Interleukin 17C ELISA Kit (orb1214197)

Description	Human Interleukin 17C, IL-17C ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	IL17C
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,765
Uniprot ID	Q9P0M4
Solubility (25°C)	9600ng/L
Sensitivity	15.49ng/L
Entrez	27189
Expiration Date	Please enquire.