

Product Datasheet

Human Fibulin 5 ELISA Kit (orb1214204)

Description	Human Fibulin 5, FBLN5 ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	FBLN5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,180
Uniprot ID	Q9UBX5
Solubility (25°C)	9600ng/L
Sensitivity	15.86ng/L
Entrez	10516
Expiration Date	Please enquire.