

Product Datasheet

Human Zinc Finger Homeobox Protein 4 ELISA Kit (orb1214205)

Description	Human Zinc Finger Homeobox Protein 4, ZFHx4 ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	ZFHx4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	393,730
Uniprot ID	Q86UP3
Solubility (25°C)	1600ng/L
Sensitivity	3.77ng/L
Entrez	79776
Expiration Date	Please enquire.