

Product Datasheet

Human Semaphorin 3A ELISA Kit (orb1214210)

Description	Human Semaphorin 3A, SEMA3A ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	SEMA3A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	88,889
Uniprot ID	Q14563
Solubility (25°C)	24ng/ml
Sensitivity	0.024ng/ml
Entrez	10371
Expiration Date	Please enquire.