

Product Datasheet

Human Histone Deacetylase 2 ELISA Kit (orb1214219)

Description	Human Histone Deacetylase 2, HDAC2 ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	HDAC2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	55,364
Uniprot ID	Q92769
Solubility (25°C)	32ng/ml
Sensitivity	0.031ng/ml
Entrez	3066
Expiration Date	Please enquire.