

Product Datasheet

Human Fibulin-1 ELISA Kit (orb1214220)

Description	Human Fibulin-1, FBLN1 ELISA Kit
Reactivity	Human
Range	3-900ug/mL
Tested Applications	ELISA
Target	FBLN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	77,214
Uniprot ID	P23142
Solubility (25°C)	960µg/mL
Sensitivity	1.54µg/mL
Entrez	2192
Expiration Date	Please enquire.