

Product Datasheet

Human Apelin ELISA Kit (orb1214222)

Description	Human Apelin, AP ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	APLN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	8,569
Uniprot ID	Q9ULZ1
Solubility (25°C)	1600ng/L
Sensitivity	3.47ng/L
Entrez	8862
Expiration Date	Please enquire.