

Product Datasheet

Human Beclin-1 ELISA Kit (orb1214223)

Description	Human Beclin-1, BECN1 ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	BECN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	51,896
Uniprot ID	Q14457
Solubility (25°C)	64ng/ml
Sensitivity	0.1ng/ml
Entrez	8678
Expiration Date	Please enquire.