

Product Datasheet

Human Ryanodine Receptor 1 ELISA Kit (orb1214227)

Description	Human Ryanodine Receptor 1, RYR1 ELISA Kit
Reactivity	Human
Range	0.1-30ng/ml
Tested Applications	ELISA
Target	RYR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	565,176
Uniprot ID	P21817
Solubility (25°C)	32ng/ml
Sensitivity	0.05ng/ml
Entrez	6261
Expiration Date	Please enquire.