

Product Datasheet

Human Interleukin 37 ELISA Kit (orb1214234)

Description	Human Interleukin 37, IL-37 ELISA Kit
Reactivity	Human
Range	7-400ng/L
Tested Applications	ELISA
Target	IL37
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,126
Uniprot ID	Q9NZH6
Solubility (25°C)	480ng/L
Sensitivity	4.56ng/L
Entrez	27178
Expiration Date	Please enquire.