

Product Datasheet

Human 1 ELISA Kit (orb1214240)

Description	Human 1, 25-Hydroxyvitamin D,1, 25-(OH)D, CYP24A1 ELISA Kit
Reactivity	Human
Range	7-1500pmol/L
Tested Applications	ELISA
Target	CYP24A1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	58,875
Uniprot ID	Q07973
Solubility (25°C)	1600pmol/L
Sensitivity	3.14pmol/L
Entrez	1591
Expiration Date	Please enquire.