

Product Datasheet

Human Podocalyxin ELISA Kit (orb1214246)

Description	Human Podocalyxin, PCX ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	PODXL
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	58,635
Uniprot ID	O00592
Solubility (25°C)	16ng/ml
Sensitivity	0.02ng/ml
Entrez	5420
Expiration Date	Please enquire.