

Product Datasheet

Human Atp-sensitive Inward Rectifier Potassium Channel 10 ELISA Kit (orb1214253)

Description	Human Atp-sensitive Inward Rectifier Potassium Channel 10, KCNJ10 ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	KCNJ10
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,508
Uniprot ID	P78508
Solubility (25°C)	4000ng/L
Sensitivity	10.32ng/L
Entrez	3766
Expiration Date	Please enquire.