

Product Datasheet

Human Progranulin ELISA Kit (orb1214259)

Description	Human Progranulin, PGRN ELISA Kit
Reactivity	Human
Range	10-700ng/ml
Tested Applications	ELISA
Target	GRN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	63,544
Uniprot ID	P28799
Solubility (25°C)	800ng/ml
Sensitivity	5.12ng/ml
Entrez	2896
Expiration Date	Please enquire.