

Product Datasheet

Human Urocortin-3 ELISA Kit (orb1214262)

Description	Human Urocortin-3, UCN3 LISA KIT
Reactivity	Human
Range	3-900ng/L
Tested Applications	ELISA
Target	UCN3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17,961
Uniprot ID	Q969E3
Solubility (25°C)	960ng/L
Sensitivity	1.56ng/L
Entrez	114131
Expiration Date	Please enquire.