

Product Datasheet

Human Chromogranin-A ELISA Kit (orb1214264)

Description	Human Chromogranin-A, CGA ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	CHGA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,688
Uniprot ID	P10645
Solubility (25°C)	9600ng/L
Sensitivity	15.21ng/L
Entrez	1113
Expiration Date	Please enquire.