

Product Datasheet

Human Caveolin-1 ELISA Kit (orb1214265)

Description	Human Caveolin-1, CAV-1 ELISA Kit
Reactivity	Human
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	CAV1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,472
Uniprot ID	Q03135
Solubility (25°C)	96ng/ml
Sensitivity	0.151ng/ml
Entrez	857
Expiration Date	Please enquire.