

Product Datasheet

Human Vimentin ELISA Kit (orb1214270)

Description	Human Vimentin, VIM ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	VIM
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	53,652
Uniprot ID	P08670
Solubility (25°C)	64ng/ml
Sensitivity	0.07ng/ml
Entrez	7431
Expiration Date	Please enquire.