

Product Datasheet

Human G1/S-specific Cyclin-D1 ELISA Kit (orb1214272)

Description	Human G1/S-specific Cyclin-D1, CYCLIN-D1 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	CCND1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	33,729
Uniprot ID	P24385
Solubility (25°C)	48ng/ml
Sensitivity	0.052ng/ml
Entrez	595
Expiration Date	Please enquire.