

Product Datasheet

Human T-lymphocyte Activation Antigen Cd86 ELISA Kit (orb1214282)

Description	Human T-lymphocyte Activation Antigen Cd86, CD86 ELISA Kit
Reactivity	Human
Range	0.1-40.ng/ml
Tested Applications	ELISA
Target	CD86
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	37,682
Uniprot ID	P42081
Solubility (25°C)	48ng/ml
Sensitivity	0.05ng/ml
Entrez	942
Expiration Date	Please enquire.