

Product Datasheet

Human Bone Morphogenetic Protein Receptor 1A ELISA Kit (orb1214285)

Description	Human Bone Morphogenetic Protein Receptor 1A, BMPR-1A ELISA Kit
Reactivity	Human
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	BMPR1A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	60,198
Uniprot ID	P36894
Solubility (25°C)	48ng/mL
Sensitivity	0.045ng/mL
Entrez	657
Expiration Date	Please enquire.