

Product Datasheet

Human Bone Morphogenetic Protein Receptor 2 ELISA Kit (orb1214286)

Description	Human Bone Morphogenetic Protein Receptor 2, BMPR2 ELISA Kit
Reactivity	Human
Range	3-900ng/L
Tested Applications	ELISA
Target	BMPR2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	115,201
Uniprot ID	Q13873
Solubility (25°C)	960ng/L
Sensitivity	1.64ng/L
Entrez	659
Expiration Date	Please enquire.