

Product Datasheet

Human Leptin ELISA Kit (orb1214289)

Description	Human Leptin, LEP ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	LEPTIN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,641
Uniprot ID	P41159
Solubility (25°C)	12.8ng/ml
Sensitivity	0.021ng/ml
Entrez	3952
Expiration Date	Please enquire.