

Product Datasheet

Human Apolipoprotein B100 ELISA Kit (orb1214291)

Description	Human Apolipoprotein B100, APO-B100 ELISA Kit
Reactivity	Human
Range	0.03-6mg/ml
Tested Applications	ELISA
Target	APOB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	515,605
Uniprot ID	P04114
Solubility (25°C)	6.4mg/ml
Sensitivity	0.01mg/ml
Entrez	338
Expiration Date	Please enquire.