

Product Datasheet

Human Semaphorin 3D ELISA Kit (orb1214297)

Description	Human Semaphorin 3D, SEMA3D ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	SEMA3D
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	89,651
Uniprot ID	095025
Solubility (25°C)	48ng/ml
Sensitivity	0.052ng/ml
Entrez	223117
Expiration Date	Please enquire.