

Product Datasheet

Human Peroxisome Proliferator-activated Receptor Gamma ELISA Kit (orb1214298)

Description	Human Peroxisome Proliferator-activated Receptor Gamma, PPARG ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	PPARG
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	57,620
Uniprot ID	P37231
Solubility (25°C)	9600ng/L
Sensitivity	14.16ng/L
Entrez	5468
Expiration Date	Please enquire.