

Product Datasheet

Human Very Low Density Lipoprotein Receptor ELISA Kit (orb1214300)

Description Human Very Low Density
Lipoprotein Receptor, VLDLR
ELISA Kit

Reactivity Human

Range 0.5-150ng/ml

**Tested
Applications** ELISA

Target VLDLR

Storage 2-8°C

Note For research use only

Application notes Sandwich

Sample Types Serum, plasma, cell culture
supernates

MW 96,098

Uniprot ID **P98155**

Solubility (25°C) 160ng/ml

Sensitivity 0.27ng/ml

Entrez **7436**

Expiration Date Please enquire.