

Product Datasheet

Human Peroxisome Proliferators Activator Receptors Alpha ELISA Kit (orb1214308)

Description	Human Peroxisome Proliferators Activator Receptors Alpha, PPAR-A ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	PPARA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	52,225
Uniprot ID	Q07869
Solubility (25°C)	9600ng/L
Sensitivity	14.22ng/L
Entrez	5465
Expiration Date	Please enquire.