

Product Datasheet

Human Dermato-Pontin ELISA Kit (orb1214309)

Description	Human Dermato-Pontin, DPT ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	DPT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,005
Uniprot ID	Q07507
Solubility (25°C)	48ng/ml
Sensitivity	0.059ng/ml
Entrez	1805
Expiration Date	Please enquire.