

Product Datasheet

Human Glypican-3 ELISA Kit (orb1214310)

Description	Human Glypican-3, GPC-3 ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	GPC3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	65,563
Uniprot ID	P51654
Solubility (25°C)	12.8ng/ml
Sensitivity	0.01ng/ml
Entrez	2719
Expiration Date	Please enquire.