

## Product Datasheet

# Human Pituitary Adenylate Cyclase Activating Polypeptide ELISA Kit (orb1214326)

|                            |                                                                           |
|----------------------------|---------------------------------------------------------------------------|
| <b>Description</b>         | Human Pituitary Adenylate Cyclase Activating Polypeptide, PACAP ELISA Kit |
| <b>Reactivity</b>          | Human                                                                     |
| <b>Range</b>               | 0.3-90ng/ml                                                               |
| <b>Tested Applications</b> | ELISA                                                                     |
| <b>Target</b>              | ADCYAP1                                                                   |
| <b>Storage</b>             | 2-8°C                                                                     |
| <b>Note</b>                | For research use only                                                     |
| <b>Application notes</b>   | Sandwich                                                                  |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates                                    |
| <b>MW</b>                  | 18,835                                                                    |
| <b>Uniprot ID</b>          | <b>P18509</b>                                                             |
| <b>Solubility (25°C)</b>   | 96ng/ml                                                                   |
| <b>Sensitivity</b>         | 0.14ng/ml                                                                 |
| <b>Entrez</b>              | <b>116</b>                                                                |
| <b>Expiration Date</b>     | Please enquire.                                                           |