

Product Datasheet

Human Calnexin ELISA Kit (orb1214332)

Description	Human Calnexin, CNX ELISA Kit
Reactivity	Human
Range	0.5-200ng/ml
Tested Applications	ELISA
Target	CANX
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	67,568
Uniprot ID	P27824
Solubility (25°C)	240ng/ml
Sensitivity	0.25ng/ml
Entrez	821
Expiration Date	Please enquire.