

Product Datasheet

Human Angiotensinogen ELISA Kit (orb1214336)

Description	Human Angiotensinogen, AGT ELISA Kit
Reactivity	Human
Range	10-1000ng/L
Tested Applications	ELISA
Target	AGT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	53,154
Uniprot ID	P01019
Solubility (25°C)	1440ng/L
Sensitivity	0.51ng/L
Entrez	183
Expiration Date	Please enquire.