

Product Datasheet

Human Clusterin ELISA Kit (orb1214342)

Description	Human Clusterin, CLU ELISA Kit
Reactivity	Human
Range	0.5-300ng/ml
Tested Applications	ELISA
Target	CLU
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	52,495
Uniprot ID	P10909
Solubility (25°C)	320ng/ml
Sensitivity	0.31ng/ml
Entrez	1191
Expiration Date	Please enquire.