

Product Datasheet

Human Aquaporin 5 ELISA Kit (orb1214369)

Description	Human Aquaporin 5, AQP-5 ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	AQP5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,292
Uniprot ID	P55064
Solubility (25°C)	16ng/ml
Sensitivity	0.021ng/ml
Entrez	362
Expiration Date	Please enquire.