

Product Datasheet

Human Oxytocin ELISA Kit (orb1214372)

Description	Human Oxytocin, OT ELISA Kit
Reactivity	Human
Range	2-600pg/ml
Tested Applications	ELISA
Target	OXT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	12,722
Uniprot ID	P01178
Solubility (25°C)	640pg/ml
Sensitivity	1.06pg/ml
Entrez	5020
Expiration Date	Please enquire.