

Product Datasheet

Human Parathyroid Hormone Related Protein ELISA Kit (orb1214376)

Description	Human Parathyroid Hormone Related Protein, PTHRP ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	PTH LH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,194
Uniprot ID	P12272
Solubility (25°C)	4000ng/L
Sensitivity	9.59ng/L
Entrez	5744
Expiration Date	Please enquire.