

Product Datasheet

Human Histone Deacetylase 8 ELISA Kit (orb1214382)

Description	Human Histone Deacetylase 8, HDAC8 ELISA Kit
Reactivity	Human
Range	0.1-35ng/mL
Tested Applications	ELISA
Target	HDAC8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	41,758
Uniprot ID	Q9BY41
Solubility (25°C)	40ng/mL
Sensitivity	0.055ng/mL
Entrez	55869
Expiration Date	Please enquire.