

Product Datasheet

Human Presenilin 2 ELISA Kit (orb1214383)

Description	Human Presenilin 2, PS-2 ELISA Kit
Reactivity	Human
Range	5-1000ng/L
Tested Applications	ELISA
Target	PSEN2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,140
Uniprot ID	P49810
Solubility (25°C)	1280ng/L
Sensitivity	2.47ng/L
Entrez	5664
Expiration Date	Please enquire.