

Product Datasheet

Human A2-Antiplasmin ELISA Kit (orb1214384)

Description	Human A2-Antiplasmin, A2-AP ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	SERPINF2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	54,566
Uniprot ID	P08697
Solubility (25°C)	2400ng/L
Sensitivity	5.07ng/L
Entrez	5345
Expiration Date	Please enquire.